

TRIANGO 60

The compact one that has what it takes



TRIANGO 60 – your new treatment & minor surgery light



TRIANGO 60 C

ceiling version



TRIANGO 60 W

wall mount



TRIANGO 60 F

roller stand



TRIANGO 60 F Battery

roller stand



TRIANGO 60	technical data
Fitted with	1x LED-Modul
Luminous power *	60 000 Lux / 1,0 m
Colour temperature *	
TRIANGO 60-1 F	fixed value 4 700 K
TRIANGO 60-3 F	changeable value 3 700 / 4 300 / 4 700 K
Lighting field 1,0 m*	ø 18 cm
Depth of illumination *	140 cm
Colour rendering index	CRI > 95
Operation equipment	Electronic converter
Voltage	100-240 V; 50/60 Hz
Luminaire body	Polyamide (PA)
Power consumption	< 30 W
Usage	Touch panel on cardanic system
Technology	Single-pole secondary
Balance of articulation joints	Spring (weight compensation)
Class of protection	I
Protection rating	IP 43 (illumination head) IP 20 (UPSmedv Pack)
Standards met	EN 60601-1, EN 60601-2-41

Special features	
TRIANGO 60-1	Continuously dimmable
TRIANGO 60-3	- Continuously dimmable 3 colour temperatures
TRIANGO Fokus 60-1	- Continuously dimmable - Focusable
TRIANGO 60-1 F Battery	- Continuously dimmable - UPSmed Pack
Mounting	- Ceiling mount - Wall mount - Roller stand
Cord length	
TRIANGO 60-1 C	Direct power supply
TRIANGO 60-1 W	Direct power supply
TRIANGO 60-1 F	3 m
TRIANGO 60-1 F Battery	3 m
Weight (net)	
TRIANGO 60-1 C	13 kg
TRIANGO 60-1 W	12,5 kg
TRIANGO 60-1 F	19 kg
TRIANGO 60-1 F Battery	31 kg
Accessories	- Sterilizable handle - Disposable cover for handle

* -10% / +20% tolerance

Versatile

versions with focusable
or colour changeable
illumination field

Homogeneous light field
thanks to unique
diamond optic

Easy cleaning

streamlined design
and the removable
handle

Hygiene made simple

good design does without
unnecessary details and
attaches importance to
closed building structures

Can be used everywhere

TRIANGO 60 is suitable
for the versatile daily
routine in family practices
and clinics as well as in
gynecology, dermatology
or veterinary medicine

**Ergonomic handling**

easy alignment for
right- and left-handed

**Intuitive operation**

thanks to an intuitive
touch panel on the
cardanic, the illumination
from 60 000 lux at 1 m
can be infinitely adjusted

Order number: TRIANGO 60-1

TRIANGO 60-1 C	D16.250.000-00813741 D16.253.000-00814312 - USA
TRIANGO Fokus 60-1 C	D16.252.000-00814311 D16.255.000-00814314 - USA
TRIANGO 60-1 W	D16.256.000-00813742 D16.274.000-00814327 - USA
TRIANGO Fokus 60-1 W	D16.258.000-00814316 D16.276.000-00814329 - USA
TRIANGO 60-1 F	D16.262.000-00813744 - CH D16.259.000-00813743 - Schuko D16.271.000-00813747 - UK D16.265.000 - USA
TRIANGO Fokus 60-1 F	D16.264.000-00814320 - CH D16.261.000-00814318 - Schuko D16.273.000-00814326 - UK D16.270.000-00814322 - USA
TRIANGO 60-1 F Battery	D16.279.000 - CH D16.278.000 - Schuko D16.282.000 - UK D16.280.000 - USA

Order number: TRIANGO 60-3

TRIANGO 60-3 C	D16.251.000-00814310 D16.254.000-00814313 USA
TRIANGO 60-3 W	D16.257.000-00814315 D16.275.000-00814328 USA
TRIANGO 60-3 F	D16.263.000-00814319 - CH D16.260.000-00814317 - Schuko D16.272.000-00814325 - UK D16.266.000-00814321 - USA





Headquarters

Derungs Licht AG
Hofmattstrasse 12
9200 Gossau
Switzerland
+41 71 388 11 66
mailbox@derungs.swiss
www.derungs.swiss



Derungs Licht AG is
a company of the
Waldmann Group

www.derungs.swiss



Register and benefit:
+ 1 year warranty and
exclusive offers

www.derungs.swiss/product-registration