





Fitted with
Illuminance field 100 cm
Illuminance 100 cm
Colour rendering index (Ra)
Work equipment
Voltage
Mains lead
Luminaire body
Weight (net)
Power consumption
Usage
Balance of articulation joints
Class of protection
Accessory
Norms
Special features

Fastening
System of protection
Medical device class

LED - Energy efficiency category A+
 ø 18 cm
 60 000 lx
 > 93
 Electronic converter
 100-240 V; 50/60 Hz
 Direct power supply
 Polyamid (PA)
 19 kg
 < 30 W
 Keypad at cardanic
 Spring
 I
 Sterilizable handle, handle cover
 EN 60601-1, EN 60601-2-41
 Stepless dimmable,
 diamond optics, sterilizable handle
 Roller stand
 IP 43 (luminaire head)
 I according to MDR

	Lighting field / 1,0 m fixed value 18 cm	lighting field / 1,0 m changeable value 18 - 25 cm	Colour temperature fixed value 4700K	Colour temperature changeable value 3700K / 4300K / 4700K	Dimmable	Order No.
TRIANGO 60-3 F CH	X			X	X	D16.263.000
TRIANGO 60-3 F Schuko	X			X	X	D16.260.000
TRIANGO 60-3 F UK	X			X	X	D16.272.000
TRIANGO 60-3 F AUS	X			X	X	D16.269.000
TRIANGO 60-3 F US	X			X	X	D16.266.000

Light distribution



Technical drawing

