




Fitted with	LED - Energy efficiency category A+ 130 000 lx
Illuminance 100 cm	> 96
Colour rendering index (Ra)	Electronic converter
Work equipment	100-240 V; 50/60 Hz
Voltage	Direct power supply
Mains lead	Polyamid (PA)
Luminaire body	15,9 kg
Weight (net)	40 W
Power consumption	Keypad at cardanic
Usage	Spring
Balance of articulation joints	I
Class of protection	Sterilizable handle, handle cover, adapter for suspended ceiling
Accessory	EN 60601-1, EN 60601-2-41
Norms	Dimmable, focusable, colour changeable, sterilizable handle
Special features	Ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm
Fastening	IP 20 according to IEC 60529 luminaire head IP 43 (horizontal position)
System of protection	I according to MDR
Medical device class	

	lighting field / 1,0 m changeable value 16 / 20 / 25 cm	Colour temperature changeable value 3700K / 4300K / 4700K	Dimmable	Order No.
TRIANGO Fokus 130-3 C	X	X	X	D16.290.000
TRIANGO Fokus 130-3 C US-HG	X	X	X	D16.297.000

Light distribution



Technical drawing

