

Treatment luminaire

TRIANGO 100-3 F D15912000-00803187

Derungs

MEDICAL LIGHTING



fitted with	LED - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 18 cm
illuminance @ 100 cm	100 000 lx
color rendering index (CRI)	> 95
colour temperatur	3700/4300/4700 K
work equipment	electronic converter
mains lead	3 m UK
luminaire body	polyamide (PA)
weight (net)	20.0 kg
Power consumption	40 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	I
Accessory	sterilizable handle
norms	EN 60601-1, EN 60601-2-41
special features	3 color shades, 4 dimming levels, Diamond optics, High-quality arm system, sterilisable grip, maintenance free
Fastening	roller stand
System of protection	IP 20 (43)
order number	D15912000 - pure white

light distribution

