

ARTro

the ART of arthroscopy



ARTro is a highly specialized electrosurgical unit dedicated for bipolar and monopolar arthroscopic procedures.

ARTro allows vaporization, cutting and bipolar coagulation in saline solution. It also offers monopolar cutting and coagulation modes which are necessary in some arthroscopic operations.

SDS system recognizes connected instruments and automatically selects the optimal working parameters. The operating modes and the output settings are selected automatically to the connected instrument.

Automatic regulation of output parameters guarantees fast and effective work.

CUT



VAPORIZATION



COAGULATION



SmartDevice System



Advantages of the device

- versatile application range: monopolar and bipolar mode
- highly simplified operation
- color display with **touchscreen**
- advanced bipolar modes for use in saline solution
- automatic adjustment of the parameters - **SpectrumResult**
- recognition of the connected devices - **SmartDevice System**
- wide range of electrodes in different shapes and sizes
- reusable accessories to reduce the costs of the procedure
- dedicated **trolley** with convenient case for accessories

SDS bipolar arthroscopic electrodes

58S-010 115 mm Bipolar needle electrode, angled 90°, 0.6 x 1.5mm, SDS cable 3m
 58S-030 170mm



58S-020 115 mm Bipolar needle electrode, angled 90°, 0.6 x 4mm, SDS cable 3m
 58S-060 170 mm



58S-040 115 mm Bipolar button electrode, angled 90°, 1.5 x 4mm, SDS cable 3m



58S-090 115 mm Bipolar loop electrode "Cobra", SDS cable 3m
 58S-011 170 mm



58S-110 115 mm Bipolar vaporization electrode, Phazer, convex, angled 70°, ball 2.4mm, SDS cable 3m
 58S-140 170 mm



58S-930 115 mm Bipolar vaporization electrode, round, VAP, SDS cable 3m
 58S-920 170 mm



Instruments with fixed cable

100-610 Electrosurgical unit ARTRO



080-061 TinyLine trolley for electrosurgical units



100-302 Two-pedal footswitch, cable 5m, 6-pin plug



See our catalogue for a comprehensive range of accessories.
 For more information, contact the regional representatives or visit www.emed.pl



Attention! This brochure does not replace the instructions for use! Refer to instructions for use!



EMED SP. Z O. O. SP. K. | Ryżowa 69a | 05-816 Opacz-Kolonia | Poland
 tel: + 48 22 723 08 00 | export@emed.pl | www.emed.pl



REF NO. 2021.06.K33.EN