

























SUPRA

SPECIÁLNÍ JEHLY

- S BROUŠENOU HRANOU JEHLY
- STERILNÍ
- NETOXICKÉ
- BEZ PYROGENU



Číslo výrobku:		Barva	Gauge	Rozměry v mm
1588 3013		žlutý	30	0,30 x 13
1588 4260		oranžový	25	0,50 x 60
1588 4450		modrý	23	0,60 x 50
1588 4460		modrý	23	0,60 x 60
1588 4470		modrý	23	0,60 x 70
1588 4480		modrý	23	0,60 x 80
1588 4580		modrý	23	0,65 x 80
1588 4650		černý	22	0,70 x 50
1588 4680		černý	22	0,70 x 80
1588 4860		zelený	21	0,80 x 60
1588 4880		zelený	21	0,80 x 80
1588 4891		zelený	21	0,80 x 100
1588 4892		zelený	21	0,80 x 120
1588 5055		žlutý	20	0,90 x 55
1588 5060		žlutý	20	0,90 x 60
1588 5080		žlutý	20	0,90 x 80
1588 5090		žlutý	20	0,90 x 90
1588 5091		žlutý	20	0,90 x 100
1588 5092		žlutý	20	0,90 x 120
1588 5260		běžový	19	1,00 x 60
1588 5291		běžový	19	1,00 x 100
1588 5292		běžový	19	1,00 x 120
1588 5492		běžový	19	1,10 x 120

Číslo výrobku:		Barva	Gauge	Rozměry v mm
1588 5650		růžový	18	1,20 x 50
1588 5660		růžový	18	1,20 x 60
1588 5675		růžový	18	1,20 x 75
1588 5680		růžový	18	1,20 x 80
1588 5691		růžový	18	1,20 x 100
1588 5930		červený	17	1,50 x 30
1588 5940		červený	17	1,50 x 40
1588 5950		červený	17	1,50 x 50
1588 5980		červený	17	1,50 x 80
1588 5990		červený	17	1,50 x 90
1588 5991		červený	17	1,50 x 100
1588 6040		šedý	16	1,60 x 40
1588 6225		šedý	15	1,80 x 25
1588 6240		šedý	15	1,80 x 40
1588 6250		šedý	15	1,80 x 50
1588 6260		šedý	15	1,80 x 60
1588 6280		šedý	15	1,80 x 80
1588 6440		zelený	14	2,00 x 40
1588 6460		zelený	14	2,00 x 60
1588 6470		zelený	14	2,00 x 70
1588 6480		zelený	14	2,00 x 80
1588 6510		zelený	14	2,00 x 100
1588 6511		zelený	14	2,00 x 120